

Small HCPs in a 12 kDa Protein Drug Analyzed by GeLC-MS/MS

Rikke Raaen Lund¹, Marie Grimstrup¹, Maheen Sayeed², Karin Abarca Heidemann², Ejvind Mørtzt¹
¹Alphalyse A/S, Odense, Denmark, ²Rockland Immunochemicals Inc., Limerick, PA, USA



Introduction

The content of low molecular weight host cell proteins (HCPs) in purified protein drugs is often difficult to evaluate, due to their low immunogenicity and poor ability to be visualized in gel-based total protein stains. The proteome of commonly used expression organisms, such as *E. coli* and Chinese Hamster cells, contains 30-40% proteins with a molecular weight below 20kDa, and these are easily missed in both in gel separations, Western blots and ELISA quantitation of the total HCP-content. To provide unbiased analysis of small as well as larger HCPs, we introduce the use of a mass spectrometry-based orthogonal method, well known from proteomics, called GeLC-MS/MS.

Here, we analyze an in process protein drug, a 12 kDa protein produced in *E. coli*, as well as the corresponding null cell lysate, using 1D-PAGE and nano-flow LC-MS/MS (GeLC-MS/MS) to achieve high coverage of small HCPs.

Figure 1: GeLC-MS/MS of null cell lysate





Figure 2: GeLC-MS/MS of the in process sample



1D-PAGE

Proteins were separated by 1D-PAGE and stained for isolation and visualization of the high concentration protein drug in separate gel fractions. Most proteins, 78% in the null cell lysate and 82% in the in process sample were only identified in a single gel fraction (light green and blue bars, fig. 1 and 2, respectively). Further, the fractionation leads to a high number of protein identifications in the fractions without the protein drug (fig 2).

Protein separation by 1D-PAGE prior to LC-MS/MS leads to identification of a high percentage of small HCPs: 46% of the HCPs identified in the null cell lysate and 37% in the in process sample were smaller than 20kDa (fig. 3, table 1 and 2). Further, the HCPs that were identified covered the entire pI range and 99% of the molecular weight range of the total *E. coli* proteome (fig. 3). Examples of small HCPs is given in table 3 and 4.

The sensitivity of the method was estimated by parallel analysis of the purified protein drug with two spiked-in standards, these standards were both identified at 50ppm.

Figure 3: Molecular weight and pI of the HCPs identified by GeLC-MS/MS and the *E. coli* proteome




Table 1: Distribution of high and low molecular weight HCPs in the null cell lysate

Protein size	HCPs identified
HMW (≥20kDa)	297
LMW (<20kDa)	256
Total	553

Table 2: Distribution of high and low molecular weight HCPs in the in process sample

Protein size	HCPs identified
HMW (≥20kDa)	96
LMW (<20kDa)	56
Total	152

Nano Flow LC-MS/MS

The proteins contained in the gel fractions were digested enzymatically into peptides using trypsin. Separation of the complex peptide mixture was achieved by nanolitre flow HPLC using a CSH (charged surface hybrid, Waters) column. This material has a high loading capacity and excellent peak shape in formic acid mobile phases. This enables highly sensitivity and accurate MS detection of low level HCPs using a qTOF mass spectrometer (Bruker Corp. fig. 4 and 5).

GeLC-MS/MS advantages

- 1D-PAGE requires very little sample preparation and has a wide mass and pI range, leading to unbiased sample analysis
- Protein digest combined with nano flow LC-MS/MS provides sensitive and accurate measurement of multiple peptides from each protein as well as sensitive and accurate measurement of peptide fragment masses

Figure 4: MS/MS spectrum of the a peptide from Ferric Uptake Regulator (in process sample)




Figure 5: CSH column and qTOF Mass spectrometer



HCP Identification Criteria

HCPs are identified by comparing the mass data to the theoretical masses of the host cell proteome. Peptides are assigned a score according to how well they match and these are summarized to a protein score. Since small proteins have fewer peptides they also have lower proteins scores even at the same molar level as larger proteins.

We have used conservative inclusion criteria that accounts for this: 1% false discovery rate, a minimum score of 25 for all peptides, and a two peptide minimum for proteins larger than 20kDa. These criteria allows highly confident HCP identification with a low false positive rate.

Table 3: Top 30 small proteins in the null cell lysate

HCP no.	Accession no.	Protein Name	Mass	pI	Score
22	tr C6E656	Ribosomal protein L9	15759	6.17	3350
29	tr C6E6C3	Ribosomal protein S7	17593	10.9	2701
31	tr C6E609	6,7-dimethyl-8-ribitylsuccinyl synthase	16147	5.15	2656
42	tr C6E600	PTS system, glucose subfamily, IIA subunit	18240	4.73	2336
46	tr C6E608	Redoxin dom-			
48	tr C6E634	S05 ribosomal			
50	tr C6E622	Inorganic diphosphate			
55	tr C6E6D0	DNA-binding			
58	tr C6E618	Histone family			
59	tr C6E6G1	HCP no.			

Table 4: Top 30 small HCPs in the in process sample

HCP no.	Accession no.	Protein Name	Mass	pI	Score		
60	tr C6E633	3	tr C6E6Q1	Ferric uptake regulator, Fur family	17012	5.68	3455
61	tr C6E603	6	tr C6E603	Peptidyle deformylase	19430	5.23	2785
62	tr C6E6E5	14	tr C6E6C1	Methionine S-sulfide reductase	15789	5.58	1881
70	tr C6E6F1	20	tr C6E634	S05 ribosomal protein L10	17757	9.04	1367
72	tr C6E632	22	tr C6E6C3	Ribosomal protein S7	17593	10.3	1082
73	tr C6E6G7	31	tr C6E6U9	Glutaredoxin	13042	4.75	906
75	tr C6E6C9	40	tr C6E6F4	Ribosomal protein L29	7269	9.38	625
76	tr C6E6C9	40	tr C6E6G3	S05 ribosomal protein S5	17534	10.23	591
79	tr C6E6W8	63	tr C6E636	S05 ribosomal protein L11	14923	9.64	399
82	tr C6E6A7	65	tr C6E6X0	DNA protection during starvation protein	18684	5.72	390
83	tr C6E6D0	70	tr C6E623	Regulator of sigma D	18288	5.65	371
86	tr C6E6D0	71	tr C6E628	Aridyldopamin synthase sulfur carrier subunit	8724	4.38	370
88	tr C6E6A6	73	tr C6E6F4	Ferritin	19468	4.77	332
96	tr C6E6A8	75	tr C6E6HB	Iron-sulfur cluster assembly scaffold protein IscU	14011	4.82	301
99	tr C6E6A12	80	tr C6E673	S05 ribosomal protein L31	8094	9.46	271
105	tr C6E6C73	81	tr C6E6E9	S05 ribosomal protein S6	15163	5.25	262
107	tr C6E6D9	83	tr C6E6B6	ATP synthase F1, delta subunit	19434	4.94	254
108	tr C6E607	88	tr C6E6F5	Uncharacterized protein GN-ECR9_4172	10323	5.18	232
110	tr C6E6G8	90	tr C6E6Q0	Flavodoxin	19896	4.21	224
112	tr C6E6F4	91	tr C6E6G2	Ribosomal protein L18	12762	10.41	222
95	tr C6E6D0	DNA-binding protein	15587	5.43	202		
96	tr C6E6G1	S05 ribosomal protein L6	18949	9.71	201		
97	tr C6E6K7	UspA domain protein	15925	6.03	200		
102	tr C6E6C0	AcyL carrier protein	8634	3.98	190		
103	tr C6E6A8	Iron-sulfur cluster assembly accessory protein	12264	4.11	189		
105	tr C6E6F5	S05 ribosomal protein L22	12219	10.23	185		
111	tr C6E6D9	Ribosomal protein S1	13950	11.39	158		
112	tr C6E6D0	PTS system, glucose subfamily, IIA subunit	18240	4.73	153		
113	tr C6E6K5	Acetyl-CoA carboxylase, biotin carboxyl carrier protein	16733	4.66	152		
115	tr C6E6K4	Thioredoxin	15887	5	147		

Conclusion

The HCPs that were identified covered 99% of the entire *E. coli* proteome in terms of molecular weight and pI. This shows that the developed GeLC-MS/MS method has no inherent limitations with respect to pI or molecular weight for HCP identification. The obtained protein identity enables an in-depth analysis of each individual HCP and a more detailed risk assessment.

HCPs Characteristics and Individual Risk Assessment

Identification of the protein names and database accession numbers enables an in-depth analysis of each individual HCP present in the drug sample. Important features for a drug risk assessment include:

- Immunological properties and presence of human B- and T-cell epitopes
- Homology to human proteins with important biological function
- Homology to the drug protein
- Enzymatic activity to modify or cleave the drug product constituents
- Hormone or hormone-like activity

Known information can be found at <http://www.uniprot.org/> including:

- Molecular function, biological process, ligand binding and cellular component
- Post translational modifications and processing
- Expression, interaction and structure
- Protein family and domains
- Sequence variations
- Publications

Prediction of protein features can be investigated at <http://www.cbs.dtu.dk/services>, including (fig. 6):

- Immunological features
- Post translational modifications
- Protein structure and function

Figure 6: Examples of services available for prediction of protein features

